Wolfram Droh GmbH

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MD

SN-0111FU 002A 23/04

Hangzhou Bever Medical Devices Co., Ltd.,
Building 2, No. 1-1, Houmuqiao,
Yongle Village, Cangqian Street,
Yuhang District,
311121 Hangzhou,
China

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(GB) Contents

1	General	Page	9
2	Safety and Responsibilities	Page	9
3	Modification of the product	Page	10
4	Before using the product	Page	10
5	Using the product	Page	11
6	Product Maintenance / Reprocessing	Page	11
7	Accessories	Page	11
8	Replacement parts and repairs	Page	11
9	Disposal	Page	11
10	Explanation of symbols	Page	11/12

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1 General

- This medical product may only be used by trained professionals under instruction of a physician.
- The product is used to create a continuous airway in the nasal area in case of obstruction or possible displacement of the nasal cavity or nasopharynx.
- The product is mainly used in emergency care.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.
- Users and/or patients within the European Union, should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Users outside of the European Union should report any serious incident that has occurred in relation to the device to the manufacturer and the regulatory authority of the country in which the user and/or patient is established.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product, please contact customer service.
- There should be sufficient experience in the use of the product in order to avoid injury to patients.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.
- The product draws air from the immediate surrounding and must therefore not be used in a toxic atmosphere.
- Use of the product in the presence of oxygen or nitrous oxide, during procedures involving lasers or electrosurgical probes, may cause the product to ignite. This can result in burns and poisoning of the patient.
- Do not apply excessive force at any time during use to avoid injury.
- A fully functioning spare of the product should always be kept ready to hand where the product is being used.
- Do not use silicone-based lubricants.

- The product must not come into contact with pointed or sharp objects to avoid cracks or other damage.
- Select the correct airway size.
- Long-term ventilation may lead to damage to the mucous membranes.
- Do not use in traumatic brain injury or midfacial fracture, skull base or ethmoid fracture, anticoagulation therapy, blood clotting disorders, transphenoidal surgery, or in fresh postoperative rhinoplasty or septoplasty.
- **Single use products must not be reused!**
They must be disposed of immediately after use in accordance with chapter 9.
 The reuse of single use products represents a potential risk for the patient or the user. It can result in contamination and/or reduced functionality.
- Contamination and/or limited functionality of the product can lead to injury, illness or death of the patient.

3 Modification of the product

- Apart from the described maintenance instructions, the product must not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- Do not expose the product to direct sunlight or other light sources.
 It must be stored and transported in a dry place and in its original packaging.
- Select the correct airway size for the patient.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ Damaged products must not be used!
 They must be disposed of in accordance with the instructions in chapter 9.
- The product contains PVC.
- The product is supplied in a sterile condition.
- To minimize contamination, gloves should be worn before and during intubation.
- The product is now ready for use.

5 Using the product

- Proceed according to the currently valid rules for emergency ventilation.
- Apply a small amount of water-soluble lubricant and with slight twisting movements gently push the tube into the nostril up to the retaining ring.
- If advancement is no longer possible, consider using the other nostril.
Do not use excessive force!
- When using a suction catheter, make sure it is the correct size and also use lubricant.
- Hold the retaining ring with one hand while carefully inserting the suction catheter with the other hand.

6 Product maintenance/Reprocessing

- Single use products must not be reused!
- After termination of the application, the product must be disposed of according to chapter 9.

7 Accessories

- The product consists of a tube and a fixed retaining ring.

8 Replacement parts and repairs

- In case of functional or material problems, contact customer service immediately.

9 Disposal

- The product must be disposed of after use.
- Please note the expiry date if applicable.
The product must not be used after the expiry date!
- Safe disposal is carried out in accordance with current national regulations.

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Single-use only / Only to be used on one patient

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